



AccuPach VI®

Item Number: 24-6200

Features:

- Table top pachymeter provides accurate corneal thickness measurements
- IOP Correction Calculation - Converts IOP measurements in seconds.
- Easy-to-use interface
- Stores information via USB memory stick
- Gentle touch probe
- Made in the USA

Sterilization Issues:

- Interchangeable Probes
- Functions with all Pach Vs and Pach VIs

Specifications:

- 65MHz Probe (sampling)
- Range of 300 to 999 microns
- Accuracy of +/-5 microns
- Resolution of +/-1 micron